



skan

Cell and Gene Therapies

Silvia Aldi, Ph.D.
Product Manager Cell and Gene Therapy
Skan AG
Silvia.aldi@skan.ch

1. Introduction to ATMPs
2. ATMP manufacturing
3. Case study I: Mesenchymal stem cell-based therapy
4. Case study II: CAR T-cell therapy
5. Fill and finish options for cell and gene therapy products
6. Summary and Conclusions

Regenerative medicine and advanced therapies financing has soared to new heights so far this year



Gene Therapy:
\$6.4B




Cell-Based IO:
\$6.6B




Cell Therapy:
\$1.1B

**\$14.1B raised in H1 2021, a 35% increase
from H1 2020**



Record Product Approvals Incoming




- Decisions are expected on 18 regenerative medicine products across 6 geographies
- Bristol Myers Squibb's Breyanzi, bluebird bio and BMS' Abecma, and Mallinckrodt's Stratagraft, all approved by the FDA, and bluebird's Skysona, approved in Europe
- Breyanzi, Abecma, and Skysona, are gene therapy/gene-modified cell therapies, which means 2021 is likely to be a record year for new approvals of this category of products
- Decisions are expected this year on four more gene therapy/gene-modified cell therapy products, with the possibility of seven total approvals which would more than double the previous record of three in 2017

Alliance for Regenerative Medicine H1 2021 report. <http://alliancerm.org/wp-content/uploads/2021/08/ARM-H1-2021-Report.pdf>


Cell & Gene Therapies: ATMPs

Advanced Therapy Medicinal Products (ATMPs) are medicines that are based on genes, tissues or cells. They represent groundbreaking new opportunities for the treatment of diseases and injuries




In vivo gene editing: direct delivery to patients:

- Viral vector
- Lipid nanoparticles
- Polymers
- Peptides



Target cells containing the transgene


Ex vivo gene editing: transgenes are delivered to cells in vitro and gene edited cells are injected into patients



Types of Therapies

- Autologous medicinal products
- Allogeneic medicinal products
- CAR T-cell therapies
- Mesenchymal stem cells-based therapy (MSCs)
- Viral vectors production
- Recombinant proteins/vaccine production
- Monoclonal antibodies/vaccine production
- CAR NK-cell therapies
- Induced pluripotent stem cells-based therapy (iPSCs)
- mRNA-based therapies (not only vaccine)
- Plasmid DNA-based therapies
- Extracellular vesicles, Exosomes

New Therapies



ATMPs developmental manufacturing process

Starting materials

- Patient/donor sample
- Viral particles
- Bacteria

Upstream process

- Cell prep
- Cell expansion
- Gene editing

Downstream

- Cell washing
- Cell separation
- Final formulation

Final Fill and finish

- Fill in final container
(glass vials, AT-closed vials, CZ vials, IV bags, bulk and nested)



ATMPs manufacturing

Starting materials

- Patient/donor sample
- Viral particles
- Bacteria

Upstream process

- Cell prep
- Cell expansion
- Gene editing

Downstream

- Cell washing
- Cell separation
- Final formulation

Final Fill and finish


- Fill in final container

Target to bring any open process into ISOLATORS




Novel therapies: manufacturing challenges

- Methodology transfer from biosafety cabinets to high aseptic work environments (isolators)
- Many different manufacturing processes/methodologies
- Wide variety of lab equipment
- Cell Therapies cannot be sterilized – high risk of contamination by open operations and growth promoting media
- Short delivery times
- Scaling-up/Scaling-out
- Stringent regulatory requirements



The complexity of Up- and Down-stream Laboratory Equipment



Customers' needs

- Aseptic controlled environment to process autologous and allogeneic therapeutics for pre-clinical and human clinical use with highest sterility assurance
- Fast transfers of biological materials in various containers (cell culture flasks, conical tubes, cryovials, bags, etc.)
- Gene manipulation by viral transduction / electroporation
- Fast and ultra-rapid decontamination cycles
- Operator and patient safety
- Acceptance from cGMP regulatory authorities (US FDA, EU EMA and other global agencies)




Examples of Cell and Gene Processes produced with isolator technology

Case Study I




Example of mesenchymal stem cell (MSCs)-based therapy



MSC

Origin: bone marrow, adipose tissue, umbelical cord

Clinical Use: for treatment of immune disorders (post-transplantation immune responses), tissue repair (damaged by stroke), cardiovascular disorder (induced ischemia), neurological and rheumatological disorders (multiple sclerosis and osteoarthritis)




<https://www.intechopen.com/chapters/69881>

Cellana Isolators for Cell and Gene Therapy


Inoculum Cell Culture

Cell thawing



Cell Harvesting

Cell harvesting and filling



Lab equipment integration and isolator design

Cellana-M solutions

- Rapid Transfer Airlock
- Customized fast decon cycles
- CO₂ Incubator
- Centrifuge
- Customized shelves




Case Study II



Case Study II


Chimeric Antigen Receptor (CAR) T-Cell Therapy: Process Overview

- **Origin:** T-cells are specific immune cells floating in the blood (cells develop in the thymus gland, T-cells)
- From donors (**allogeneic**) or patients (**autologous**)
- **Biological specifications:** Chimeric Antigen Receptor (CAR) T-cells are genetically manipulated *in vitro* to target specific cancer antigens expressed on cancer cells to treat lymphomas, leukemia




<https://www.nature.com/articles/s41434-021-00246-w.pdf>

CAR T-cell Manufacturing Process in Isolator



Cellana-L

- High Flexibility
- Interchangeable modules
- Scalability and Versatility



Automated or Manual ATMPs filling

Key features

- Filling technology minimizing the contamination risks
- Containers must be resistant to cryopreservation agent (DMSO)
- For safe cryogenic storage at -80°C and in the nitrogen liquid tank



- Ready-to-fill closed vial
- Scalable filling equipment
- Full validation package
- Particularly suitable for Cell and Gene Therapies



AT-Closed Vial[®]
Technology



<https://www.aseptictech.com/>

Crystal PURE M1



Cell and Gene Isolators

Cellana Crystal L-1



IV Bag Filling-small scale

Summary and Conclusions

- As lab-based processes are introduced for clinical and GMP production, increased controls are required
- Isolator technology is a natural next step to increase process and quality control
- Challenges in ergonomics and transfers must be addressed, including deco. practices and equipment integration
- End goal is not necessarily perfection, but improved quality



skan

Thank you for your time!

Silvia Aldi, PhD
Product Manager Cell and Gene Therapy
silvia.aldi@skan.ch

