



A FORUM FOR SUPPORTERS OF THE DISTRIBUTED NETWORK PROTOCOL

Distributed Network Protocol

This certificate confirms that the device described below has successfully complied with the requirements as outlined by the DNP 3.0 Technical Committee's DNP 3.0 Conformance Testing procedure. This certificate entitles the bearer to promote the below stated device as having fully passed the DNP 3.0 Subset Definitions v2.6 rev1.

Manufacturer: Siemens


Device Model: SIPROTEC 4, 7SC80

Firmware Number: V01.00.05

Level Successfully Passed: Level 2

Testing House: 

Date: November 23, 2012



Tester:

A handwritten signature in blue ink, appearing to read 'P.H.S. Ermens'.

P.H.S. Ermens