

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

**Siemens AG, IC SG EA PRO D 2 4
Wernerwerkdamm 5, 13629 Berlin, Germany**

the Certificate No: **Z10482** for the PROFINET IO Device:

Model Name: SIPROTEC 4
Revision: SW/FW: V01.02.05; HW: 3
GSD: GSDML-V2.25-Siemens-SIPROTEC4-20140307.xml
Identnumber: 0x002A; 0x0E03
DAP: IDD_1: SIPROTEC4 with EN100-E, 0x00000001

This certificate confirms that the product has successfully passed the certification tests with the following PROFINET scope:

<input checked="" type="checkbox"/> Hardware	Auto_Negotiation, Auto_Polarity, Auto_Crossover
<input checked="" type="checkbox"/> Conformance Class B	RT_CLASS_1, RTA, LLDP, SNMP, MIB-II, LLDP-MIB

Test Report Number: PN186-2
Authorized Test Laboratory: SIEMENS AG, Fürth, Germany


The tests were executed in accordance with the following documents:
 "Test Specifications for PROFINET IO devices, Version 2.2.4 from December 2010",
 "Test Cases for PN-Tester for PROFINET IO devices, Version 2.2.4.16.21".

This certificate is granted according to the document:
 "Framework for testing and certification of PROFIBUS and PROFINET products".
 For all products that are placed in circulation by May 02, 2017 the certificate is valid for life.


Board of PROFIBUS Nutzerorganisation e. V.



(Official in Charge)



(Karsten Schneider)

(K.-P. Lindner)