



GBS INVENTORY RECOMMENDER

Half empty? Half full? Well balanced.

There are some things that better don't go out of stock. On the other hand, a warehouse with crammed shelves isn't what you're looking for either. What both have in common: They cost a lot of money that you could invest in your business. The good news is, you can handle your inventory intelligently from now on with the GBS Inventory Recommender.

The GBS Inventory Recommender is a cloud-based tool that taps directly into your Enterprise Resource Planning (ERP) system and uses artificial intelligence to transform inventory-related data into knowledge that helps you getting closer to just-in-time manufacturing.

With this smart inventory management solution, you can increase the accuracy of lead times by up to 53%, reduce delayed deliveries by up to 35% and lower stock levels to achieve up to 50% inventory related savings.* At the same time, you can increase visibility over supplier performance and process bottlenecks.

Striking the balance between keeping your supply chain lean, always having just enough stock in your warehouse and still being able to react fast and flexible to your clients' requirements? Possible.

With the GBS Inventory Recommender, you're one decisive step closer to turn this ambition into reality.

Intelligent inventory management – made possible by Siemens Global Business Services.

*Source: GBS Inventory Recommender pilot project

Siemens AG

Global Business Services
Otto-Hahn-Ring 6
81739 München

[siemens.com/gbs](https://www.siemens.com/gbs)

SIEMENS